


GOAT ANTI-GSTA3 ANTIBODY

SKU: EB13063



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	glutathione S-transferase A3 glutathione S-aryltransferase A3 glutathione S-alkyltransferase A3 glutathione S-alkyltransferase A3 S-(hydroxyalkyl)glutathione lyase A3 GST class-alpha member 3 GTA3 GSTA3-3 glutathione S-transferase alpha 3 GSTA3
Accession ID	NP_000838.3
Blocking Peptide	EBP13063
Immunogen	Peptide with sequence C-RKPPADAKALE, from the C Terminus of the protein sequence according to NP_000838.3.
Peptide Sequence	C-RKPPADAKALE
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	2940
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 24kDa band observed in Human Placenta lysates (calculated MW of 25.3kDa according to NP_000838.3). Recommended concentration: 0.1-0.3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

