

GOAT ANTI-GZMA ANTIBODY

SKU: EB09312



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	granzyme A Hanukah factor serine protease Granzyme A (Cytotoxic T-lymphocyte-associated serine esterase-3 CTL tryptase HFSP CTLA3 granzyme A (granzyme 1, cytotoxic T-lymphocyte-associated serine esterase 3) GZMA
Accession ID	NP_006135.1
Blocking Peptide	EBP09312
Immunogen	Peptide with sequence C-QLTEKAKINKY, from the internal region of the protein sequence according to NP_006135.1.
Peptide Sequence	C-QLTEKAKINKY
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	3001
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)