



Telephone: (650) 697-3600



GOAT ANTI-GZMA ANTIBODY

SKU: EB09312

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / granzyme A|Hanukah factor serine protease)|Granzyme A (Cytotoxic T-lymphocyte-associated serine

Alias esterase-3|CTL tryptase|HFSP|CTLA3|granzyme A (granzyme 1, cytotoxic T-lymphocyte-associated serine esterase

Names 3)|GZMA

Accession

NP 006135.1 ID

Blocking

EBP09312 **Peptide**

Immunogen NP_006135.1. Peptide with sequence C-QLTEKAKINKY, from the internal region of the protein sequence according to

Peptide

C-QLTEKAKINKY Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human

Human

3001 Gene ID

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

ELISA

Detection

Antibody detection limit dilution 1:8000.

Limit

Application Pep-ELISA **Type**