



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-HAX1 ANTIBODY**

**SKU:** EB09160



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

250kDa 150kDa 100kDa 75kDa 50kDa

25kDa 20kDa

15kDa





Telephone: (650) 697-3600

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

HAX1|HCLS1 associated protein X-1|FLJ17042|FLJ18492|FLJ93803|HCLSBP1|HS1BP1|SCN3|HCLS1 (and PKD2) **Alias** 

Names

associated protein|HS1 binding protein|OTTHUMP00000034190

**Accession** 

NP\_006109.2; NP\_001018238.1

**Blocking** 

EBP09160 **Peptide** 

Peptide with sequence C-TRHEADSSPRGDPES, from the internal region of the protein sequence according to **Immunogen** 

NP\_006109.2; NP\_001018238.1.

**Product Comments** 

This antibody is expected to recognize both reported isoforms (NP 006109.2 and NP 001018238.1).

**Peptide** 

Sequence

**C-TRHEADSSPRGDPES** 

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

**Predicted** 

**Species** 

Human

Reactive

Human **Species** 

Human

10456

Gene ID

**Product** 

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

**ELISA** 

**Detection** 

Antibody detection limit dilution 1:2000.

Limit

Western Approx 33kDa band observed in Human Kidney, Liver and Testis lysates (calculated MW of 31.6kDa according

to NP\_006109.2). Recommended concentration: 1-3µg/ml. Blot

**Application Type** 

Pep-ELISA, WB

## **GALLERY IMAGES**



Telephone: (650) 697-3600



250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa 20kDa

15kDa