

# GOAT ANTI-HAX1 ANTIBODY

**SKU:** EB09160

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

## SPECIFICATIONS

<b>Formulation</b>	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
<b>Unit Size</b>	100 µg
<b>Storage Instructions</b>	Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Synonym / Alias Names</b>	HAX1 HCLS1 associated protein X-1 FLJ17042 FLJ18492 FLJ93803 HCLSBP1 HS1BP1 SCN3 HCLS1 (and PKD2) associated protein HS1 binding protein OTTHUMP00000034190
<b>Accession ID</b>	NP_006109.2; NP_001018238.1
<b>Blocking Peptide</b>	EBP09160
<b>Immunogen</b>	Peptide with sequence C-TRHEADSSPRGDPEs, from the internal region of the protein sequence according to NP_006109.2; NP_001018238.1.
<b>Product Comments</b>	This antibody is expected to recognize both reported isoforms (NP_006109.2 and NP_001018238.1).
<b>Peptide Sequence</b>	C-TRHEADSSPRGDPEs
<b>Purification Method</b>	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
<b>Shipping Instructions</b>	Refrigerated
<b>Predicted Species</b>	Human
<b>Reactive Species</b>	Human
<b>Human Gene ID</b>	10456
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a>
<b>ELISA Detection Limit</b>	Antibody detection limit dilution 1:2000.
<b>Western Blot</b>	Approx 33kDa band observed in Human Kidney, Liver and Testis lysates (calculated MW of 31.6kDa according to NP_006109.2). Recommended concentration: 1-3µg/ml.
<b>Application Type</b>	Pep-ELISA, WB

## DOCUMENTS

- [Data Sheet](#)

## GALLERY IMAGES

