

GOAT ANTI-HCFC1 ANTIBODY

SKU: EB10460



SPECIFICATIONS

Unit Size	100 µg
Synonym /	
Alias	HCFC1 VCAF MGC70925 host cell factor C1 (VP16-accessory protein) host cell factor 1 HFC1 HCF-1 HCF1 CFF
Names	
Accession ID	NP_005325.2
Blocking Peptide	EBP10460
Immunogen	Peptide with sequence CGPSPSCLVQ, from the internal region of the protein sequence according to NP_005325.2.
Peptide Sequence	CGPSPSCLVQ
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Human Gene ID	3054
Mouse Gene ID	15161
Rat Gene ID	363519
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Preliminary experiments gave an approx 150kDa band in lysates of Human Peripheral Blood Mononucleocytes and of cell line Jurkat after 1µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated size of 209kDa according to NP_005325.2. The 150kDa band was successfully blocked by incubation with the immunizing peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates? Have any further splice variants/modified forms been reported?
Application Type	Pep-ELISA