

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-HDGFRP3 (AA89-103) ANTIBODY

SKU: EB12826



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

HDGFRP3|HDGF-2|HDGF2|HRP-3|hepatoma-derived growth factor 2|hepatoma-derived growth factor-related

protein 3 Names

Accession

NP_057157.1

Blocking

Peptide

EBP12826

Immunogen

Peptide with sequence C-ENNPGVKFTGYQAIQ, from the internal region of the protein sequence according to

NP 057157.1.

Peptide

Sequence

C-ENNPGVKFTGYQAIQ

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Method using the immunizing peptide.

Shipping

Instructions

Refrigerated

Human

50810

Predicted

Species

Human, Dog, Cow

Reactive

Species

Human

Gene ID

Product

Grade

ELISA

Detection

Limit Western

Antibody detection limit dilution 1:128000.

Blot

Approx 25kDa band observed in Human Brain (Amygdala) lysates (calculated MW of 22.6kDa according to NP_057157.1). Recommended concentration: 0.3-1µg/ml.

Application

Type

Pep-ELISA, WB

GALLERY IMAGES



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa