

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-HE4 (RAT) ANTIBODY

SKU: EB12490

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Re4|Wfdc2 Alias

Names

Accession NP 775132.1 ID

Blocking

EBP12490 **Peptide**

Peptide with sequence C-RTSPLSRGQVSTKPP, from the internal region of the protein sequence according to Immunogen

NP 775132.1.

Peptide

C-RTSPLSRGQVSTKPP Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Rat **Species**

Rat Gene ID 286888

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:128000.

Limit

Western **Blot**

Preliminary experiments in Rat Colon, Heart, Kidney and Uterus lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any

results been reported with other antibodies/lysates?

Application

Pep-ELISA **Type**