

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-HER6 (ZEBRAFISH) ANTIBODY**

**SKU:** EB08009

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias her6|hairy-related 6|HER-6|cb129

**Names** 

Accession

NP 571154.1

**Blocking** 

EBP08009 **Peptide** 

Peptide with sequence KRKIKQTVTSLPHQ-C, from the N Terminus of the protein sequence according to Immunogen

NP 571154.1.

**Peptide** 

KRKIKQTVTSLPHQ-C Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

**Species** 

Zebrafish

**Product** 

Grade

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. a mazon aws.com/assets/products/image/aspiring\_medium.png$ 

**ELISA** 

**Detection** Limit

Antibody detection limit dilution 1:128000.

Western

Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the product is investigative grade. We

Blot would appreciate any feedback from people in the field - have any results been reported with other

antibodies/lysates?

**Application Type** 

Pep-ELISA