

GOAT ANTI-HES7 ANTIBODY

SKU: EB10783



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	transcription factor HES-7 mHes7 bHLH factor Hes7 OTTMUSP00000006176 bHLHb37 hairy and enhancer of split 7 (Drosophila) Hes7
Accession ID	NP_149030.2
Blocking Peptide	EBP10783
Immunogen	Peptide with sequence PRSPGQDAEALASCY, from the internal region of the protein sequence according to NP_149030.2.
Peptide Sequence	PRSPGQDAEALASCY
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Pig, Cow
Human Gene ID	84667
Mouse Gene ID	84653
Rat Gene ID	305811
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)