

GOAT ANTI-HIP1R ANTIBODY

SKU: EB05496



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	huntingtin interacting protein-1-related KIAA0655 HIP1R HIP3 huntingtin interacting protein 12 HIP12
Names	
Accession ID	NP_003950.1
Blocking Peptide	EBP05496
Immunogen	Peptide with sequence C-DKKDGIYPAQLVNY, from the C Terminus of the protein sequence according to NP_003950.1.
Peptide Sequence	C-DKKDGIYPAQLVNY
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat
Reactive Species	Human
Human Gene ID	9026
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 120-125kDa band seen in Human Liver lysate (Expected MW from amino acid sequence according to XP_208650 is 116kDa). Recommended for use at 2-4µg/ml.