

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-HIP1R/ HIP12 ANTIBODY**

**SKU:** EB09186

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

huntingtin interacting protein-1-related|huntingtin interacting protein **Alias** 

**Names** 

12|MGC47513|KIAA0655|ILWEQ|HIP3|HIP12|FL|14000|huntingtin interacting protein 1 related|HIP1R

Accession

NP 003950.1 ID

**Blocking** 

EBP09186 **Peptide** 

Peptide with sequence C-RQDHQLDKKDGIYP, from the C Terminus of the protein sequence according to Immunogen

NP 003950.1.

**Peptide** Sequence

C-RQDHQLDKKDGIYP

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

**Species** 

Human, Dog, Cow, Pig

Human

9026 Gene ID

**Product** Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

**ELISA** 

**Detection** 

Antibody detection limit dilution 1:32000.

Limit

**Application** Pep-ELISA **Type**