

GOAT ANTI-HIP1R/ HIP12 ANTIBODY

SKU: EB09186



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	huntingtin interacting protein-1-related huntingtin interacting protein 12 MGC47513 KIAA0655 ILWEQ HIP3 HIP12 FLJ14000 huntingtin interacting protein 1 related HIP1R
Accession ID	NP_003950.1
Blocking Peptide	EBP09186
Immunogen	Peptide with sequence C-RQDHQLDKKDG IYP, from the C Terminus of the protein sequence according to NP_003950.1.
Peptide Sequence	C-RQDHQLDKKDG IYP
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Dog, Cow, Pig
Human Gene ID	9026
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)