



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-HIP2 ANTIBODY

SKU: EB05750



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa

100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa



Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

UBE2K|HIP2|huntingtin interacting protein 2|LIG|HYPG|ubiquitin-protein ligase|ubiquitin carrier Synonym / protein|ubiquitin-conjugating enzyme E2-25 KDA|ubiquitin-conjugating enzyme E2K (UBC1 homolog,

yeast)|DKFZp564C1216|DKFZp686J24237|E2-25K Names

Accession

Alias

ID

NP_005330.1; NP_001104582.1; NP_001104583.1

Blocking

EBP05750 **Peptide**

Peptide with sequence C-SWDVETATELLLSN, from the C Terminus of the protein sequence according to **Immunogen**

NP_005330.1; NP_001104582.1; NP_001104583.1.

Product This antibody is expected to recognise isoform 1 (NP_005330.1), isoform 2 (NP_001104582.1) and isoform 3

Comments (NP_001104583.1).

Peptide

C-SWDVETATELLLSN Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse, Dog, Pig, Cow **Species**

Reactive

Human **Species**

Human

3093

Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade

In paraffin embedded Human Brain shows pixulate cytoplasm staining in neurons of the Cortex. Recommended **IHC Results**

concentration, 1-2µg/ml.

ELISA

Antibody detection limit dilution 1:2000. Detection

Limit

Western Approx 24kDa band observed in Jurkat cell lysates (calculated MW of 22.4 kDa according to NP_005330).

Blot Recommended for use at 1-3 µg/ml.

Application

Pep-ELISA, WB, IHC **Type**



Telephone: (650) 697-3600

GALLERY IMAGES

250kDa 150kDa

100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

