

GOAT ANTI-HMHA1 / HA-1 ANTIBODY

SKU: EB05126



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	histocompatibility (minor) HA-1 HMHA1 KIAA0223 HLA-HA1 minor histocompatibility antigen HA-1 HA-1
Accession ID	NP_036424.2
Blocking Peptide	EBP05126
Immunogen	Peptide with sequence TLGSCRERQPEFV, from the C Terminus of the protein sequence according to NP_036424.2.
Product Comments	To avoid confusion please note that this antibody is designed to recognise the whole protein as coded by NP_036424.2 and won't recognise the peptide found eluted from HLA-A*0201 molecules by den Haan et al (Science. 1998 Feb 13;279(5353):1054-7, PMID: 9461441)
Peptide Sequence	TLGSCRERQPEFV
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	23526
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)