

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-HMHA1 / HA-1 ANTIBODY

SKU: EB05126

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias histocompatibility (minor) HA-1|HMHA1|KIAA0223|HLA-HA1|minor histocompatibility antigen HA-1|HA-1

Names

Accession NP 036424.2 ID

Blocking

EBP05126 **Peptide**

Immunogen Peptide with sequence TLGSCRERQPEFV, from the C Terminus of the protein sequence according to NP 036424.2.

To avoid confusion please note that this antibody is designed to recognise the whole protein as coded by **Product**

NP_036424.2 and won't recognise the peptide found eluted from HLA-A*0201 molecules by den Haan et al Comments

(Science. 1998 Feb 13;279(5353):1054-7, PMID: 9461441)

Peptide

TLGSCRERQPEFV Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human Species

Human

23526

Gene ID Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:32000.

Limit

Application Pep-ELISA **Type**