


GOAT ANTI-HOXA10 ANTIBODY

SKU: EB12204



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names homeobox protein Hox-A10|homeobox protein Hox-1H|homeobox protein Hox-1.8|homeobox protein HOXA10|homeobox protein 1H|homeo box A10|PL|HOX1H|HOX1.8|HOX1|homeobox A10|HOXA10

Accession ID NP_061824.3

Blocking Peptide EBP12204

Immunogen Peptide with sequence PQATSCSFAQNIKE, from the internal region of the protein sequence according to NP_061824.3.

Peptide Sequence PQATSCSFAQNIKE

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species	Human, Mouse, Pig, Cow
Reactive Species	Human
Human Gene ID	3206
Mouse Gene ID	15395
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Testis. Recommended concentration: 7µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:2000.
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

