




Telephone: (650) 697-3600

GOAT ANTI-HOXA10 ANTIBODY

SKU: EB12204

e erest



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym / homeobox protein Hox-A10|homeobox protein Hox-1H|homeobox protein Hox-1.8|homeobox protein **Alias** HOXA10|homeobox protein 1H|homeo box A10|PL|HOX1H|HOX1.8|HOX1|homeobox A10|HOXA10 **Names**

Accession

NP_061824.3

Blocking

ID

EBP12204 **Peptide**

Immunogen

Peptide with sequence PQATSCSFAQNIKE, from the internal region of the protein sequence according to NP_061824.3.

Peptide

PQATSCSFAQNIKE

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

e erest

Predicted

Human, Mouse, Pig, Cow

Species Reactive

Human **Species**

Human

3206 **Gene ID**

Mouse **Gene ID** 15395

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. a mazon aws.com/assets/products/image/elite_medium.png$ Grade

IHC Results Paraffin embedded Human Testis. Recommended concentration: 7µg/ml.

ELISA

Antibody detection limit dilution 1:2000. Detection

Limit

Application Pep-ELISA, IHC Type

GALLERY IMAGES

