


GOAT ANTI-HOXA5 (AA83-92) ANTIBODY

SKU: EB12199



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names homeobox protein Hox-A5|homeobox protein Hox-1C|homeobox protein HOXA5|homeo box A5|homeo box 1C|HOX1C|HOX1.3|HOX1|homeobox A5|HOXA5

Accession ID NP_061975.2

Blocking Peptide EBP12199

Immunogen Peptide with sequence C-EPRYSQPATS, from the internal region of the protein sequence according to NP_061975.2.

Peptide Sequence C-EPRYSQPATS

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human, Mouse, Pig, Cow

Reactive Species Human, Mouse, Rat, Pig

Human Gene ID 3202

Mouse Gene ID 15402

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA Detection Limit Antibody detection limit dilution 1:8000.

Western Blot Approx 28kDa band observed in Human, Mouse, Rat and Pig Kidney lysates (calculated MW of 29.2kDa according to Mouse NP_034583.1). Recommended concentration: 1-3µg/ml.

Application Type Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

