

GOAT ANTI-HOXB4 ANTIBODY

SKU: EB09005



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	homeo box B4 homeo box 2F HOX2F HOX2 HOX-2.6 homeobox B4 HOXB4
Accession ID	NP_076920.1
Blocking Peptide	EBP09005
Immunogen	Peptide with sequence C-EYSQSDYLPSDH, from the internal region of the protein sequence according to NP_076920.1.
Peptide Sequence	C-EYSQSDYLPSDH
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Human Gene ID	3214
Mouse Gene ID	15412
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)