

GOAT ANTI-HOXD12 ANTIBODY

SKU: EB11908



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	homeobox protein Hox-D12 homeobox protein Hox-4H homeo box D12 Hox-4.7, mouse, homolog of HOX4H homeobox D12 HOXD12
Accession ID	NP_067016.3
Blocking Peptide	EBP11908
Immunogen	Peptide with sequence C-PTAKDGPEEQAK, from the internal region of the protein sequence according to NP_067016.3.
Peptide Sequence	C-PTAKDGPEEQAK
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	3238
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:1000.
Western Blot	Approx 37kDa band observed in Human Heart lysates (calculated MW of 29.0kDa according to NP_067016.3). This molecular weight is routinely observed by other sources. Recommended concentration: 2-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

