Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-HRASLS / H-REV107 (INTERNAL) ANTIBODY

SKU: EB08325

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias H-REV107 protein-related protein|HSD28|HRASLS1|H-REV107|A-C1|HRAS-like suppressor|HRASLS

Names

Accession NP 065119.1 ID

Blocking

EBP08325 Peptide

Peptide with sequence C-RYGEGVSEQANR, from the internal region of the protein sequence according to Immunogen

NP 065119.1.

Peptide

C-RYGEGVSEQANR Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse, Rat **Species**

Human

57110

Gene ID Mouse

Gene ID

27281

Rat Gene ID 288025

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

ELISA

Detection Antibody detection limit dilution 1:2000.

Limit

Western Preliminary experiments in Human Heart, Skeletal Muscle, and Testis lysates gave no specific signal but low

Blot background (at antibody concentration up to 1µg/ml).

Application

Pep-ELISA Type





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600