

GOAT ANTI-HRASLS / H-REV107 (INTERNAL) ANTIBODY

SKU: EB08325



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	H-REV107 protein-related protein HSD28 HRASLS1 H-REV107 A-C1 HRAS-like suppressor HRASLS
Accession ID	NP_065119.1
Blocking Peptide	EBP08325
Immunogen	Peptide with sequence C-RYGEGVSEQANR, from the internal region of the protein sequence according to NP_065119.1.
Peptide Sequence	C-RYGEGVSEQANR
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat
Human Gene ID	57110
Mouse Gene ID	27281
Rat Gene ID	288025
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Western Blot	Preliminary experiments in Human Heart, Skeletal Muscle, and Testis lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml).
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)