


GOAT ANTI-HRASLS3 / HREV107-3 ANTIBODY

SKU: EB08338



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names adipose-specific PLA2|HRAS-like suppressor 3|MGC118754|HREV107-3|HREV107|HRASLS3|HREV107-1|AdPLA|phospholipase A2, group XVI|PLA2G16

Accession ID NP_009000.2

Blocking Peptide EBP08338

Immunogen Peptide with sequence RAPIPEPKPGDLIE-C, from the N Terminus of the protein sequence according to NP_009000.2.

Product Comments Reported variants represent identical protein:NP_009000.2 and NP_001121675.1).

Peptide Sequence RAPIPEPKPGDLIE-C

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Cow
Reactive Species	Human
Human Gene ID	11145
Mouse Gene ID	225845
Rat Gene ID	24913
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Brain (Cortex). Recommended concentration: 5µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:128000.
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

