


GOAT ANTI-HTR2A ANTIBODY

SKU: EB08522



SPECIFICATIONS

Unit Size	100 µg
Synonym /	
Alias	5-HT2 receptor HTR2 5-HT2A RP11-147L20.2 5-hydroxytryptamine (serotonin) receptor 2A
Names	
Usage	Flow Cytometry: Flow cytometric analysis of A549 cells. Recommended concentration:
Summary	10ug/ml.
Accession	
ID	NP_000612.1
Blocking	
Peptide	EBP08522
Immunogen	Peptide with sequence KNSKQDAKTTDNDCS, from the internal region of the protein sequence according to NP_000612.1.
Peptide	
Sequence	KNSKQDAKTTDNDCS

Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	3356
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:1000.
Application Type	Pep-ELISA, FC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

