
Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-HTR2A ANTIBODY**





## **SPECIFICATIONS**

100 μg **Unit Size** 

Synonym /

5-HT2 receptor|HTR2|5-HT2A|RP11-147L20.2|5-hydroxytryptamine (serotonin) receptor 2A Alias

**Names** 

<strong>Flow Cytometry:</strong> Flow cytometric analysis of A549 cells. Recommended concentration: Usage

**Summary** 10ug/ml.

Accession

NP\_000612.1 ID

**Blocking Peptide** 

EBP08522

Peptide with sequence KNSKQDAKTTDNDCS, from the internal region of the protein sequence according to Immunogen NP\_000612.1.

**Peptide** Sequence

KNSKQDAKTTDNDCS





Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

Shipping Refrigerated Instructions

**Predicted** 

Human **Species** 

Reactive

Human

**Species** Human

3356

Gene ID **Product** 

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. a mazon aws.com/assets/products/image/elite\_medium.png$ 

Grade **ELISA** 

Detection

Antibody detection limit dilution 1:1000.

Limit

Type

 $\begin{array}{c} \textbf{Application} \\ \textbf{Pep-ELISA, FC} \end{array}$ 

## **GALLERY IMAGES**

