

GOAT ANTI-HUMAN IGG ANTIBODY HRP CONJUGATED

SKU: EB2ND-50-3001-HRP



SPECIFICATIONS

Formulation	Phosphate buffer saline, pH 7.2, with 10 mg/ml BSA and 0.01% Gentamicin sulphate as preservative. Quantity: 0.5ml.
Unit Size	500 µg
Storage Instructions	Store at 4°C for short durations. Store at -20°C for longer periods (> 1 month). Minimize freeze thaw cycles.
Usage Summary	Suitable for immunoblot (western or dotblot), ELISA, and other peroxidase-based immunoassays. ELISA: 1:50000 WB: 1:25000-1:50000. Please optimize dilutions for your particular assay conditions. Starting test dilution: 0.5µg/ml for ELISA and WB.
Immunogen	Purified whole human IgG (Horseradish Peroxidase/ HRP conjugated).
Purification Method	Ammonium sulphate precipitation of serum followed by antigen affinity purification.
Shipping Instructions	Refrigerated
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png