




Telephone: (650) 697-3600

GOAT ANTI-HUMAN C-REACTIVE PROTEIN ANTIBODY (HRP) ANTIBODY

SKU: EB13101



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Synonym /

Alias PTX1|CRP|C-reactive protein

Names

Accession

NP 000558.2

Peptide with sequence Purified Human C reactive protein, from the of the protein sequence according to **Immunogen**

NP 000558.2.

Peptide

Purified Human C reactive protein Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping Instructions

Refrigerated

Predicted

Species

Human

Reactive

Human **Species**

Human

1401 Gene ID

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection Antibody detection limit dilution 1:128000.

Limit

Western Approx. 25kDa band observed in Human Liver and Serum Lysates. Recommended concentration: 0.3µg/ml.

Blot Primary incubation was 1 hour.

Application

Pep-ELISA, WB **Type**

GALLERY IMAGES

250kDa 150kDa

100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600