


GOAT ANTI-HUMAN C-REACTIVE PROTEIN ANTIBODY (HRP) ANTIBODY

SKU: EB13101



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias PTX1|CRP|C-reactive protein

Names

Accession ID NP_000558.2

Immunogen Peptide with sequence Purified Human C reactive protein, from the of the protein sequence according to NP_000558.2.

Peptide Sequence Purified Human C reactive protein

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human

Reactive Species Human

Human Gene ID 1401

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection Limit Antibody detection limit dilution 1:128000.

Western Blot Approx. 25kDa band observed in Human Liver and Serum Lysates. Recommended concentration: 0.3µg/ml. Primary incubation was 1 hour.

Application Type Pep-ELISA, WB

GALLERY IMAGES

