

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-HUMAN IGG ANTIBODY HRP CONJUGATED

SKU: EB2ND-50-3001-HRP

SPECIFICATIONS

Formulation 0.5 ml. Phosphate buffer saline, pH 7.2, with 10 mg/ml BSA and 0.01% Gentamicin sulphate as preservative. Quantity:

500 μg **Unit Size**

Instructions Store at 4°C for short durations. Store at -20°C for longer periods (> 1 month). Minimize freeze thaw cycles.

Suitable for immunoblot (western or dotblot), ELISA, and other peroxidase-based immunoassays. ELISA:

Usage Summary

1:50000 WB: 1:25000-1:50000. Please optimize dilutions for your particular assay conditions. Starting test

dilution: 0.5µg/ml for ELISA and WB.

Immunogen Purified whole human IgG (Horseradish Peroxidase/ HRP conjugated).

Purification Method

Ammonium sulphate precipitation of serum followed by antigen affinity purification.

Shipping

Instructions

Refrigerated

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png