

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-PERILIPIN 1 (C TERMINUS) ANTIBODY

SKU: EB13103



Telephone: (650) 697-3600



250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

lipid droplet-associated protein|PERI|perilipin|Plin|Perilipin1|Perilipin 1|PLIN1

Names

Accession

NP_002657.3

Blocking

Peptide

EBP13103

Peptide with sequence PAPRLAQPRRSC, from the C Terminus of the protein sequence according to **Immunogen**

NP 002657.3.

Product

This is an alternative product to EB07728. **Comments**

Peptide

Sequence

PAPRLAQPRRSC

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Instructions

Refrigerated

Human

Predicted

Species

Reactive

Human **Species**

Human

5346

Gene ID **Product**

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite medium.png

ELISA

Detection

Western

Antibody detection limit dilution 1:2000.

Limit

Blot

Type

Approx 55-58kDa band observed in Human Adipose tissue lysates (calculated MW of 56.0kDa according to Human NP_002657.3). Recommended concentration: 0.03-0.1ug/ml. Primary incubation 1 hour at room

temperature.

Application

Pep-ELISA, WB

GALLERY IMAGES



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa