

GOAT ANTI-PERILIPIN 1 (C TERMINUS) ANTIBODY

SKU: EB13103



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	lipid droplet-associated protein PERI perilipin Plin Perilipin1 Perilipin 1 PLIN1
Accession ID	NP_002657.3
Blocking Peptide	EBP13103
Immunogen	Peptide with sequence PAPRLAQPRRSC, from the C Terminus of the protein sequence according to NP_002657.3.
Product Comments	This is an alternative product to EB07728.
Peptide Sequence	PAPRLAQPRRSC
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	5346
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Western Blot	Approx 55-58kDa band observed in Human Adipose tissue lysates (calculated MW of 56.0kDa according to Human NP_002657.3). Recommended concentration: 0.03-0.1ug/ml. Primary incubation 1 hour at room temperature.
Application Type	Pep-ELISA, WB

GALLERY IMAGES

