

GOAT ANTI-PERILIPIN 1 (C TERMINUS) ANTIBODY

SKU: EB13102

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias Perilipin 1|Perilipin1|Peri|PERIA|Plin|FPLD4|PLIN1|PLIN|PERI

Names

Accession ID NP_002657.3

Blocking Peptide EBP13102

Immunogen Peptide with sequence RVSDSFFRPSVMEPC, from the C Terminus of the protein sequence according to NP_002657.3.

Product Comments This is an alternative product to EB07728.

Peptide Sequence RVSDSFFRPSVMEPC

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human

Reactive Species Human

Human Gene ID 5346

Mouse Gene ID 103968

Rat Gene ID 25629

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection Limit Antibody detection limit dilution 1:128000.

Western Blot Approx 55-58kDa band observed in Human Adipose tissue lysates (calculated MW of 56.0kDa according to Human NP_002657.3). Recommended concentration: 0.1-0.3µg/ml. Primary incubation 1 hour at room temperature.

Application Type Pep-ELISA, WB

GALLERY IMAGES

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

