

GOAT ANTI-HYALURONIDASE 1 (MOUSE, RAT) ANTIBODY

SKU: EB07872



SPECIFICATIONS

Unit Size	100 µg
Synonym /	
Alias	hyaluronidase 1 Hyal-1 Hyal hyaluronoglucosaminidase 1
Names	
Accession ID	NP_032343.2
Blocking Peptide	EBP07872
Immunogen	Peptide with sequence SDKTSTKESCQ, from the internal region of the protein sequence according to NP_032343.2.
Peptide Sequence	SDKTSTKESCQ
Shipping Instructions	Refrigerated
Predicted Species	Mouse, Rat
Mouse Gene ID	15586
Rat Gene ID	367166
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:1000.
Western Blot	Preliminary experiments in Mouse and Rat Brain lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA