


GOAT ANTI-ICAM1 ANTIBODY

SKU: EB11753



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names P3.58|major group rhinovirus receptor|intercellular adhesion molecule 1 (CD54), human rhinovirus receptor|intercellular adhesion molecule 1|ICAM-1|ICAM1|human rhinovirus receptor|cell surface glycoprotein P3.58|CD54|BB2

Usage Summary Immunofluorescence: Strong expression of the protein seen in the plasma membrane and cytoplasm of U251 cells. Recommended concentration: 10µg/ml.

Accession ID NP_000192.2

Blocking Peptide EBP11753

Immunogen Peptide with sequence SNCPDGQSTAKT, from the internal region of the protein sequence according to NP_000192.2.

Product	The immunizing peptide represents part of the extracellular domain, and it does not overlap any known
Comments	glycosite.
Peptide Sequence	SNCPDGQSTAKT
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	3383
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:4000.
Application Type	Pep-ELISA, IF

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

