


Telephone: (650) 697-3600

GOAT ANTI-ICAM1 ANTIBODY

SKU: EB11753



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

100 μg **Unit Size**

Storage Instructions

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / P3.58|major group rhinovirus receptor|intercellular adhesion molecule 1 (CD54), human rhinovirus

Alias

receptor|intercellular adhesion molecule 1|ICAM-1|ICAM1|human rhinovirus receptor|cell surface glycoprotein

Names P3.58|CD54|BB2

Usage Summary Immunofluorescence: Strong expression of the protein seen in the plasma membrane and

cytoplasm of U251 cells. Recommended concentration: 10µg/ml.

Accession

NP_000192.2

Blocking Peptide

EBP11753


Immunogen

Peptide with sequence SNCPDGQSTAKT, from the internal region of the protein sequence according to

NP 000192.2.











Product The immunizing peptide represents part of the extracellular domain, and it does not overlap any known

Comments glycosite.

Peptide

SNCPDGQSTAKT

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method

using the immunizing peptide.

Shipping Instructions

Refrigerated

Predicted

Species

Human

Reactive **Species**

Human

Human

3383 Gene ID

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection

Antibody detection limit dilution 1:4000.

Limit

Application

Type

Pep-ELISA, IF

GALLERY IMAGES

