

GOAT ANTI-ICBP90 / UHRF1 ANTIBODY

SKU: EB05943



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	UHRF1 ubiquitin-like, containing PHD and RING finger domains, 1 Np95 ICBP90 FLJ21925 transcription factor ICBP90 RNF106 MGC138707 RING finger protein 106 inverted CCAAT box-binding protein of 90 kDa nuclear zinc finger protein Np95 transcription factor ICBP90 ubiquitin-like with PHD and ring finger domains 1 hNP95 E3 ubiquitin-protein ligase UHRF1
Accession ID	NP_001041666.1; NP_037414.3
Blocking Peptide	EBP05943
Immunogen	Peptide with sequence C-LNQLFPGYGNGR, from the C Terminus of the protein sequence according to NP_001041666.1; NP_037414.3.
Product Comments	This antibody is expected to recognise both reported isoforms.
Peptide Sequence	C-LNQLFPGYGNGR
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog
Reactive Species	Human
Human Gene ID	29128
Mouse Gene ID	18140
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Western Blot	Approx 100kDa band observed in nuclear lysates of cell line Jurkat (calculated MW of 91.1kDa according to NP_037414.3). The observed molecular weight corresponds to earlier findings in literature with different antibodies (Hopfner, Cancer Res. 2000 Jan 1;60(1):121-8. PMID: 10646863). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

