

GOAT ANTI-ICEBERG ANTIBODY

SKU: EB05229



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	ICEBERG caspase-1 inhibitor ICEBERG
Accession ID	NP_067546
Immunogen	Peptide with sequence C-EEDPQLASKMGLH, from the C Terminus of the protein sequence according to NP_067546.
Peptide Sequence	C-EEDPQLASKMGLH
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	59082
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:1000.
Western Blot	Western Blot: Band observed in Mouse Heart extracts at approx 25kDa with 1 hour incubation of primary.