

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-ICEBERG ANTIBODY**

**SKU:** EB05229

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias ICEBERG caspase-1 inhibitor|ICEBERG

Names

Accession NP 067546

Immunogen Peptide with sequence C-EEDPQLASKMGLH, from the C Terminus of the protein sequence according to NP\_067546.

**Peptide** C-EEDPQLASKMGLH Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

**Shipping** 

Instructions Refrigerated

**Predicted** 

Human **Species** 

Human

59082 **Gene ID** 

**Product** https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png

Grade **ELISA** 

**Detection** Antibody detection limit dilution 1:1000.

Limit

Western

Western Blot: Band observed in Mouse Heart extracts at approx 25kDa with 1 hour incubation of primary. Blot