



Telephone: (650) 697-3600



GOAT ANTI-IL-10 (AA25-39) ANTIBODY

SKU: EB11757

SPECIFICATIONS

Unit Size

100 μg

Synonym / Alias

TGIF|T-cell growth inhibitory factor|OTTHUMP00000034567|MGC126451|MGC126450|interleukin-10|interleukin

10|IL10A|IL-10|IL10|cytokine synthesis inhibitory factor|CSIF Names

Accession

ID

NP_000563.1

Blocking

Peptide

EBP11757

Immunogen

Peptide with sequence QSENSCTHFPGNLPN, from the internal region of the protein sequence according to

NP 000563.1.

Peptide

Sequence

QSENSCTHFPGNLPN

Shipping

Instructions

Refrigerated

Human

3586

Predicted

Species

Human

Gene ID

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

ELISA

Detection Limit

Antibody detection limit dilution 1:32000.

Western

Blot

Preliminary experiments in all Human cell line lysates tested gave no specific signal but low background at

antibody concentration up to 1µg/ml.

Application

Type

Pep-ELISA