

GOAT ANTI-IL-10 (AA25-39) ANTIBODY

SKU: EB11757



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	TGIF T-cell growth inhibitory factor OTTHUMP00000034567 MGC126451 MGC126450 interleukin-10 interleukin 10 IL10A IL-10 IL10 cytokine synthesis inhibitory factor CSIF
Accession ID	NP_000563.1
Blocking Peptide	EBP11757
Immunogen	Peptide with sequence QSENSCTHFPGNLPN, from the internal region of the protein sequence according to NP_000563.1.
Peptide Sequence	QSENSCTHFPGNLPN
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	3586
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Preliminary experiments in all Human cell line lysates tested gave no specific signal but low background at antibody concentration up to 1µg/ml.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)

