

GOAT ANTI-IL-15 (AA70-82) ANTIBODY

SKU: EB11758

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Instructions	
Synonym / Alias Names	OTTHUMP00000228805 OTTHUMP00000228804 OTTHUMP00000218093 OTTHUMP00000164617 OTTHUMP00000164616 MGC9721 interleukin-15 interleukin 15 IL-15 IL15
Accession ID	NP_000576.1; NP_751915.1
Blocking Peptide	EBP11758
Immunogen	Peptide with sequence C-DATLYTESDVHPS, from the internal region of the protein sequence according to NP_000576.1; NP_751915.1.
Product Comments	This antibody is expected to recognize both reported isoforms (NP_000576.1; NP_751915.1).
Peptide Sequence	C-DATLYTESDVHPS
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	3600
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 23kDa band observed in Human Bone Marrow lysates (calculated MW of 18.1kDa according to NP_000576.1). The observed molecular weight corresponds to the glycosylated form. Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

