

GOAT ANTI-IL-17A (AA126-137) ANTIBODY

SKU: EB11728

250kDa
150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	OTTHUMP00000016597 interleukin 17A interleukin 17 (cytotoxic T-lymp IL-17A IL17A IL-17 IL17 cytotoxic T-lymphocyte-associated serine esterase 8 cytotoxic T-lymphocyte-associated protein 8 cytotoxic T-lymphocyte-associated antigen 8 CTLA-8 CTLA8
Accession ID	NP_002181.1
Blocking Peptide	EBP11728
Immunogen	Peptide with sequence PPHCPNSFRLEK, from the internal region of the protein sequence according to NP_002181.1.
Peptide Sequence	PPHCPNSFRLEK
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	3605
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 19kDa band observed in Human Lymph Nodes and Thymus lysates (calculated MW of 17.5kDa according to NP_002181.1). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES



250kDa
150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa