

GOAT ANTI-IL-17D (AA55-65) ANTIBODY

SKU: EB12085

250kDa
150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias interleukin-27|interleukin-17D|interleukin 27|IL27|IL-27|IL-22|IL-17D|interleukin 17D|IL17D

Names

Accession ID NP_612141.1

Blocking Peptide EBP12085

Immunogen Peptide with sequence C-HHTLQLGPREQ, from the internal region of the protein sequence according to NP_612141.1.

Peptide Sequence C-HHTLQLGPREQ

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human, Mouse

Reactive Species Human

Human Gene ID 53342

Mouse Gene ID 239114

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA Detection Limit Antibody detection limit dilution 1:16000.

Western Blot Approx 20 and 23kDa bands observed in Human Brain (Frontal Cortex) lysates (calculated MW of 21.9kDa according to NP_612141.1). The observed 23kDa band corresponds to the glycosylated form. Recommended concentration: 1-3µg/ml.

Application Type Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

