


GOAT ANTI-IL-25 / IL-17E (AA66-79) ANTIBODY

SKU: EB12023



250kDa
150kDa
100kDa
75kDa

50kDa

37kDa

25kDa
20kDa

15kDa

SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	IL25 interleukin 25 IL17E interleukin 17E interleukin-17E interleukin-25
Names	
Accession ID	NP_073626.1; NP_758525.1
Blocking Peptide	EBP12023
Immunogen	Peptide with sequence RPNRHPESCRASED, from the internal region of the protein sequence according to NP_073626.1; NP_758525.1.
Product Comments	This antibody is expected to recognize both reported isoforms (NP_073626.1; NP_758525.1).
Peptide Sequence	RPNRHPESCRASED
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	64806
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 23+19kDa bands observed in Human Prostate lysates (calculated MW of 20.3kDa according to NP_073626.1 and of 18.5kDa according to NP_758525.). Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

