

GOAT ANTI-IL-25 / IL-17E ANTIBODY

SKU: EB12024

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	IL25 interleukin 25 IL17E interleukin 17E interleukin-17E interleukin-25
Accession ID	NP_073626.1; NP_758525.1
Blocking Peptide	EBP12024
Immunogen	Peptide with sequence HGEKGTHKGYCLER, from the internal region of the protein sequence according to NP_073626.1; NP_758525.1.
Product Comments	This antibody is expected to recognize both reported isoforms (NP_073626.1; NP_758525.1).
Peptide Sequence	HGEKGTHKGYCLER
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	64806
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 21kDa band observed in Human Prostate lysates (calculated MW of 20.3kDa according to NP_073626.1). Recommended concentration: 0.5-1.5µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES



