


GOAT ANTI-IL-6 (AA138-148) ANTIBODY

SKU: EB11755



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias OTTHUMP00000198490|OTTHUMP00000198489|OTTHUMP00000198486|OTTHUMP00000158544|interleukin-6|interleukin BSF-2|interleukin 6 (interferon, beta 2)|interferon beta-2|IL-6|IL6|IFN-beta-2|IFNB2|hybridoma growth factor|HSF|HGF|CTL differentiation factor|CDF|BSF-2|BSF2|B-cell stimulatory factor 2|B-cell differentiation factor|B cell stimulatory factor-2

Usage Summary Immunofluorescence: Strong expression of the protein seen in U2OS cells. Recommended concentration: 10µg/ml.

Accession ID NP_000591.1; NP_001358025.1; NP_001305024.1

Blocking Peptide EBP11755

Immunogen Peptide with sequence C-EQARAVQMSTK, from the internal region of the protein sequence according to NP_000591.1; NP_001358025.1; NP_001305024.1.

Peptide Sequence C-EQARAVQMSTK

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	3569
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Tonsil. Recommended concentration: 5µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:64000.
Application Type	Pep-ELISA, IF, IHC

GALLERY IMAGES

