




GOAT ANTI-IL12A ANTIBODY

SKU: EB09028



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	natural killer cell stimulatory factor 1, 35 kD subunit interleukin-12 alpha chain interleukin 12, p35 cytotoxic lymphocyte maturation factor 1, p35 NF cell stimulatory factor chain 1 IL-12, subunit p35 P35 NKSF1 NFSK IL-12A CLMF interleukin 12A (natural killer cell stimulatory factor 1, cytotoxic lymphocyte maturation factor 1, p35) IL12A
Accession ID	NP_000873.2
Blocking Peptide	EBP09028
Immunogen	Peptide with sequence PCTSEEDHEDITK, from the internal region of the protein sequence according to NP_000873.2.
Peptide Sequence	PCTSEEDHEDITK
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated



Predicted Species	Human, Dog, Pig
Reactive Species	Human
Human Gene ID	3592
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	In paraffin embedded Small Intestine shows granular staining of paneth cells. Recommended concentration: 4-6µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:16000.
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

