

GOAT ANTI-IL17C (AA160-172) ANTIBODY

SKU: EB12090

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	IL17C interleukin 17C CX2 IL-17C IL-21 cytokine CX2 interleukin-17C
Accession ID	NP_037410.1
Blocking Peptide	EBP12090
Immunogen	Peptide with sequence RRPCSRDGSGGLPT, from the internal region of the protein sequence according to NP_037410.1.
Peptide Sequence	RRPCSRDGSGGLPT
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	27189
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 20kDa band observed in Human Lymph Node lysates (calculated MW of 21.8kDa according to NP_037410.1). Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa