

GOAT ANTI-IL33 (AA165-178) ANTIBODY

SKU: EB12001

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	IL33 interleukin 33 C9orf26 DVS27 IL1F11 NF-HEV NFEHEV RP11-575C20.2 DVS27-related protein IL-1F11 IL-33 interleukin-1 family member 11 interleukin-33 nuclear factor for high endothelial
Names	venules nuclear factor from high endothelial venules
Accession ID	NP_254274.1; NP_001186569.1; NP_001186570.1
Blocking Peptide	EBP12001
Immunogen	Peptide with sequence C-ESQHPSNESGDGVD, from the internal region of the protein sequence according to NP_254274.1; NP_001186569.1; NP_001186570.1.
Product Comments	This antibody is expected to recognize all reported isoforms (NP_254274.1; NP_001186569.1; NP_001186570.1).
Peptide Sequence	C-ESQHPSNESGDGVD
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	90865
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:4000.
Western Blot	Approx 30kDa band observed in lysates of cell line Jurkat (calculated MW of 30.7kDa according to NP_254274.1). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

